



Economical Freeze Dryer



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific

Darsun Scientific






Economical Freeze Dryer

DS3L-60A-2 DS3L-60B-2 DS3L-60C-2
DS3L-60E-2 DS3L-60D-2



There is another appearance Pioneer 3L, with name Pioneer -3L -60-2, same specification only different appearance.



Model					
	DS3L-60A-2	DS3L-60B-2	DS3L-60C-2	DS3L-60D-2	DS3L-60E-2
condenser temperature	-80°C				
Ice condenser capacity	3L/24h				
Dimension (mm)	393X479X451(without chamber)				
Vacuum pressure (pa)	< 1Pa (standard < 5Pa)				
vacuum pump	2XZ-2(2L/s)				
*optional vacuum pump	DRV-10(10 ³ /h)				
Refrige system	2 compressors Cascade				
Defrost	Heating coil				
Rated power	1350w				
Net weight (Appro)	45kg(without vacuum pump, vacuum pump 20kg)				

