



# Manual Sealing Machine

---



## India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
[info@darsunscientific.in](mailto:info@darsunscientific.in)  
(+91) 9999136670, 7835864003

## USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific



Darsun Scientific

# Manual Sealing Machine

DS100-S



Model	DS100-S	
Technical Parameters	1.	Sealing Width: Max. 300mm ;
	2.	Sealing Line Width: 10mm ;
	3.	Temperature: 60°C~200°C adjustable;
	4.	AC Power: 110V /220V±10% 50/60Hz
	5.	Maximum Current: 3.2A;
	6.	Power: 500W
	7.	Humidity: ≤90% (R.H); Atmospheric Pressure: 50KPa~106KPa;
	8.	Size: 450×390×240 (L×W×H) mm
	9.	Net Weight: 9Kgs
	10.	The device is used only by turning the handle, without additional pressure.
	11.	It is strictly prohibited to seal Tyvek bags (120°C low temperature bags)!
	12.	Plug: Optional for European / American
	13.	Manual sealing, manual cutting, not recommended for disinfection supply room (not applicable to Tyvek sterilization bag), suitable for: laboratories, private clinics, etc.
	14.	Alloy heating components, temperature adjustable, ultra-high temperature automatic protection, low temperature automatic alarm, overcurrent protection and other leading design with fast heating speed, stable performance, beautiful appearance, safe usage, easy operation, low repair rate advantages. Built-in double-edge cutting knife to ensure smooth incision.
	15.	DS100-S manual sealing machine is a special equipment for paper-plastic bags, paper-plastic three-dimensional bags and paper-paper bags that need to be heat sealed on site. It can cut and seal paper plastic bags below 300mm, which mainly meets the sealing of paper plastic bags before sterilization of various medical items in hospitals, medical consumables production units or pharmaceutical factories, and can meet the needs of high temperature steam sterilization, ethylene oxide and gamma-ray sterilization.