



Stackable incubator shaker (Button control)



India Office
16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters
30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States

DarsunScientific



Darsun Scientific

Stackable incubator shaker (Button control)

DS-S3/DS-EDS-4 series




Stackable incubator shaker are the idea choice for optimizing laboratory space without compromising the advantages of performing large capacity incubation the product are designed with sliding platform are light controller for easy of access the samples

Model	DS-S3	DS-S3L	DS-S4
Description	Stackable Incubator Shaker		
Speed	0rpm, 30-300rpm		
Speed Accuracy	1rpm		
Orbit Diameter	Ø26mm		
Max. Capacity	Clamp Capacity:-50mlx84 100mlx42/ 200mlx40/ 500mlx28/ 1000mlx15/ 2000mlx12/3000mlx8/5000mlx6 Sticky mat :-50mlx124 100mlx93 150x68 200mlx54 250mlx52 300mlx40 500mlx28 1000mlx15 2000mlx12 3000mlx8 5000mlx6	Clamp Capacity:-520mlx84/100mlx42 /150mlx42 200mlx40 /250mlx40/500mlx28/ 1000mlx15 /2000mlx12 3000mlx8 Sticky mat:- 50mlx124 100mlx93 150mlx68 200mlx54 250mlx52 300mlx40 500mlx28 1000mlx15 2000mlx12 3000mlx8	Clamp Capacity:- 520mlx84/100mlx42 /150mlx42 200mlx40 /250mlx40/500mlx28/ 1000mlx15 /2000mlx12 3000mlx8 Sticky mat:- 50mlx124 100mlx93 150mlx68 200mlx54 250mlx52 300mlx40 500mlx28 1000mlx15 2000mlx12
Platform Dimensions	740×500mm		
Temp. Range	15°C Below Ambient-60°C (Min. 4°C)		
Temp. Accuracy	0.1°C		
Temp. Uniformity	±0.3°C (@37°C)	±0.8°C (@37°C)	±0.3°C (@37°C)
Light Strength	~20000Lux (up to 10cm from platform)		
Auto-Defrost	Yes		
Auto-run after Power Recovery	Yes		
Display	LCD		
Timer	0-999h59min		
Outside Dimensions	1240×777×670mm		1240×777×570mm
Operation Chamber Dimension	400mm	363mm	300mm
Power	<1100W	<1300W	<1100W
Power Supply	AC220V±22V, 50 / 60Hz		
Net Weight	220kg (Single Unit)		200kg (Single Unit)
Maximum Stackable Height	2m (3 Units)		1.7m (3 Units)



Category	Details
Part Number	DS-S3IS-S4
Platform Dimensions	740x500 mm
140x140 mm Sticky Mat Quantity	15 pcs



Part number	Specifications
DS-A13	136x95mm



96 well microplate holder

Part number	Specifications	Note
DS-A29	984x630x225 mm	Used for DS-S3 DS-S4 SERIES



Rack for stacking incubator shakers

Part number	Specifications	Note
DS-A116	864x500x43mm	IS-4 series



Static Shelf

Part number	Specifications	Note
DS-A81	652x232mm	For DS-4 series



Window shade

Part number	Specifications	Note
DS-A107	130x98x85	For 24 well plate Dimension :- 127x85x51 mm (With cover and sticky mat) Used for DS-S3, DS-4 Series



24-Deep Well Plate Clamp

Part number	Specifications	Note
DS-A33	500ml	Made of stainless steel can be fixed of incubator shaker



Part number	Specifications	Note
DS-A34	1000ml	Made of stainless steel can be fixed of incubator shaker

