



WB

# Water Bath



#### India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

#### USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific



Darsun Scientific

# Water Bath

DS-1-WB100/DS-2-WB100  
 DS-4-WB100/DS-6-WB100  
 DS-8F-WB100



## FEATURES

 One-piece stainless steel tank

 High corrosion resistance

 High precise temperature control

 Digital display thermostatic

Model	DS-1-WB100	DS-2-WB100	DS-4-WB100	DS-6-WB100	DS-8F-WB100
Temperature Range	Room temp-100°C				
Power	300W	600W	800W	1200W	1500W
Accuracy	±0.1°C				
Working Chamber Size (mm)	160*160*130	300*150*140	300*300*140	450*300*130	620*300*120