



**WB**

## **Water Bath With Magnetic Stirrer**

---



### **India Office**

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

### **USA Headquarters**

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States



Darsun Scientific

## Water Bath With Magnetic Stirrer

DS-WB-6S/DS-WBS-6P



### FEATURES



Brushless motor



Strong Stirring Power



PID temperature control



Stainless steel inner tank

Model	DS-WB-6S	DS-WBS-6P
Speed / Increment	/	100-1200RPM / 50RPM
Temperature Range	Room Temp. +5°C -100°C	Room Temp. +5°C -100°C
Liner Size	200*200*130mm(depth)-2L flask or 3L beaker	
Product Size (mm)	260*265*325	260*265*325
Display Accuracy	0.1°C	0.1°C
Power	800W	
Timing Range	1min-23h59min	