



Automatic Cutting Machine



India Office

16 S/F B/S, T-Block Extn.
Jain Colony Part-2 Uttam Nagar,
New Delhi-110059
info@darsunscientific.in
(+91) 9999136670, 7835864003

USA Headquarters

30 N Gould St, Ste R
Sheridan, Wyoming 82801
United States


DarsunScientific



Darsun Scientific

Automatic Cutting Machine

DS200-S/DS200-SL



Equipment Model	DSMY200-S		DS200-SL
Control System	4.3" color LCD touch screen		
Cutting Width	$\leq 350\text{mm}$ (Multiple rolls can be cut simultaneously)		$\leq 520\text{mm}$ (multiple rolls can be cut simultaneously)
Cutting Length	$\geq 50\text{mm}$		
Cutting Speed	$18 \pm 0.5\text{m/min}$		
AC Power	$110\text{V}/220\text{V} \pm 10\% 50/60\text{Hz}$		
Counting Function	Can count 0–9999		
Cutting Quantity	Can cut 1–999		
Cutting Mode	Left and right circulation cutting, more beautiful and neat		
Blade	Round blade, can achieve multi-point use		
Optional Accessories	Extended roller-able workbench, sealing machine workstation		
Application	Suitable for paper-plastic sterilization bags, Tyvek sterilization bags, and other roll bags		
Tailing Monitoring	When only the last part of the bag is left, the equipment will automatically monitor and suspend to prevent the tailing from entering the equipment		
Function Type	Automatic paper feeding and cutting; set required length and quantity for automated completion. Includes counting function and a floating elastic pressure rolling structure for various bag thicknesses		
Multiple Rolls Cutting	Multiple rolls of sterilization bags with different widths can be loaded at one time for simultaneous cutting		
Parameter Setting	Built-in clock, parameter settings, and automatic storage function		
Self-check Function	Automatic detection during working process and automatic alarm prompt in case of failure		