



# Continuous Medical Sealing Machine

---



## India Office

16 S/F B/S, T-Block Extn.  
Jain Colony Part-2 Uttam Nagar,  
New Delhi-110059  
info@darsunscientific.in  
(+91) 9999136670, 7835864003

## USA Headquarters

30 N Gould St, Ste R  
Sheridan, Wyoming 82801  
United States

**Darsun**Scientific



Darsun Scientific

# Continuous Medical Sealing Machine

DSY100-A/DSY100-AT



Model	DSY100-A	DSY100-AT
Control System	LED digital display	4.3" color LCD touch screen
Sealing Margin	0-35mm adjustable	
Sealing Speed	10 ±0.5min (can be updated to 18m/min)	
AC Power	110V/220V\10%,50/60Hz	
Sealing Line Width	12mm	
Working Temperature	60-220°C set arbitrarily	
Temperature Deviation	≤1%	
Function Type	Continuous sealing, sealing counting	
Heat and Cool	Room temperature rises to 180°C in less than 40s, auxiliary cooling design	
Parameter Setting	Built-in clock, parameter settings, automatic storage function	
Safety Guarantee	The deviation between the sealing temperature and the set working temperature exceeds 4°C, it will automatically stop working	
Heating Method	Advanced flat heating element, can dry burn, high temperature resistance, long service life, high thermal efficiency	
Housing Material	Carbon steel / stainless steel / aluminium alloy	
Indentation Style	Flat wheel sealing (Default) / Single indentation sealing / Pattern wheel sealing	
Optional Accessories	Extended roller-able workbench, medical sealing machine workstation	
Application	Paper-plastic sterilization bags, Tyvek sterilization bags (1073B, 1059B), paper-plastic three-dimensional bags, PE bags, aluminum foil bags, etc.	